# DEI® Clever fibre COMPOSITE

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[www.deitalia.it](http://www.deitalia.it)

## DEI Clever Fibre COMPOSITE

HIGH RESISTANCE UNIVERSAL PHOTO-POLYMERIZING HYBRID COMPOSITE

Product of dental use. Keep away from children. Date of the specifications: 11/2017

**Available packages**  
DEI Clever Fibre Composite: One 4 gr syringe  
- LIGHT COMFC 2  
- MEDIUM COMFC 3

### Characteristics and classification

DEI Clever Fibre Composite is a product planned and developed for offering exceptional performances of resistance in restorations of areas subjected to strong masticatory stresses. Compared to values shown by other composites, presence of a special filler, **glass microfibres**, particular sinterization techniques used in phase

of production, high filler percentage (73% in weight) give DEI Clever Fibre Composite resistance three times higher.

Also in relation to polymerization shows value of great importance (1% in volume). Being to these properties, the technical performances of DEI Clever Fibre Composite are not obtainable at present by other marketed composites.

DEI Clever Fibre Composite is available in two colorings: light and medium.

DEI Clever Fibre Composite meets the ISO 4049 regulations.

### Indications

- 1° and 2° class restorations
- Restorations of cavities in sandwich technique
- Prosthetic structures for bridges and crowns

### Composition

DEI Clever Fibre Composite is a light-cured dental restorative material, based on dimethacrylate resins, containing inorganic fillers in different sizes, fillers, stabilizers and the photoinitiator system (CQ, DMAEMA)

responsive to visible blue light (400-500 nm).

### Way of utilization

#### Preparation of the cavity

First of all, treat the enamel with water and pumice; isolate with dam the operating area. The preparation must be not contaminated by saliva or other materials. Prepare the cavity seeking for retaining as much enamel as possible. To improve the closing of the edges and the final safety of the restoration, chamfer at 45°-60° the edges of the preparation. Preparing rear teeth, avoid to leave enamel without support of dentin.

In deep cavities we recommend the application of Calcium Hydroxide (not included in the package) for protection of the pulp.

### Technique of etching and adhesive preparation

See DEI New Extra Bond user instructions.

### Application of composite

Apply layers of composite not thicker than 3 mm. Light cure for 40 seconds.

Repeat the procedure up to full restoration of the tooth.

### Finishing technique

Before the photopolymerization of last layer of composite, apply DEI Surface Hardener and cure for 40 seconds, wash and dry. The hardening of the composite shall be complete. Using this product avoids the necessity of eliminating the superficial layer of inhibition with burs or other products, with shorter work times.

Should increases of appearance and resistance of the construction be wanted, apply DEI Seal Coat Fast. Cure for 60 seconds; let polymerize for other 20 seconds. Take away the superficial layer of DEI Seal Coat Fast using abrasive paste and small brush made with goat hair, small brush of cotton and DEI Polish Paste. After finishing, an extremely hard, resistant and shiny surface is obtained. Further, this technique considerably decreases the build-up of tartar and bacterial plaque on the restored tooth.

### Cautions

Restorations with composite don't use products containing eugenol. Keep refrigerated at constant temperature of 4°C or in environment with fresh and constant temperature lower than 25°C. Take the pack from the refrigerator some hours before use. Don't expose to sources of light or heat.

After use seal again the syringe with its suitable cap. To release the pressure of the material, turn the piston halfway antitockwise. Avoid contact with the skin, as in case of contact with the skin, use the specific product. For the utilization of the accessory products, read always the instructions for use of the specific product.

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