DEI® NEW EXTRA BOND TWO

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AVAILABLE PACKAGES:

Introducing Kit (REF. NEB K2) contains:

- One 5 ml bottle of **DEI® New Extra Primer**
- One 5 ml bottle of **DEI® New Extra Bond Two**
- One 5 ml bottle of **DEI® New Extra Catalyst**
- Two 2 ml syringes of **DEI® Etching Gel**
- Working accessories

REFILL PACKAGES:

- One 5 ml bottle of **DEI® New Extra Primer** (REF. NEB P)
- One 5 ml bottle of **DEI® New Extra Bond Two** (REF. NEB 2)
- One 5 ml bottle of **DEI® New Extra Catalyst** (REF. NEB C)
- Four 2 ml syringes of **DEI® Etching Gel** (REF. MORD 3)

INDICATIONS

1. Direct and indirect restorations, cementing and core build-ups with light-curing composites based on methacrylates.
2. Intraoral restorations of ceramic veneers, fillings in composite, ceramic restorations without additional primer.
3. Treatment of cervical dentinal hypersensitivity and/or radicular post-scaling and root planing.
4. Protective liner concrete glass ionomer fillings.
5. Sealing of cavities before a restoration with amalgam.
6. Direct and indirect restorations with self-curing or dual-curing composites based on methacrylates (in combination with DEI® New Extra Catalyst, activator for the dual-curing).
7. Cementing and core build-ups with all self-curing or dual-curing composites.

8. Cementation of root posts with self-curing or dual-curing composites.
9. Cementation of zirconia veneers with all self-curing or dual-curing composites.
10. Cementing crowns and bridges, metal-free ceramic and zirconia.

ADVANTAGES

- Free choice of the type of etching.
 - 1) Self-etch or thanks to the formulation of DEI® New Extra Bond Two it is not necessary to etch dentin as the product, without the use of added acid, promotes adhesion autonomously.
 - 2) Selective-enamel-etch.
 - 3) Total-etch-and-rinse.
- Free choice of the type of polymerization (with light or dual).
- Chemical and mechanical adhesion.
- Bond values and high longevity.
- Resistance to hydrolytic degradation.
- Reduces the risk of post-operative sensitivity.
- Flexible and secure. Suitable for every clinical situation.
- Suitable for direct and indirect restorations.
- Adhesion to various materials, such as metal, ceramic, disilicate lithium without additional primers for direct interventions on the patient.
- Multi component for a reliable adhesion and superior durability.

CHARACTERISTICS

DEI® New Extra Bond Two is a multi-component adhesive with selective polymerization usable in combination with all the techniques of etching: self-etch, selective-enamel-etch and total-etch-and-rinse. It is also indicated for both restorations direct and indirect, and for cementing with all composites light-curing or dual-curing.

DEI® New Extra Bond Two comes up in in three parts (primer, adhesive and catalyst). DEI® New Extra Bond Two can be directly applied on dentin without the need of phase separate etching phase with phosphoric acid, while getting the same adhesive features of the total-etch etchings. DEI® New Extra Bond Two tolerates the presence of residual moisture, it can also be used in a universal way with any etching technique both in enamel that in dentin (self-etch, selective-enamel-etch and total-etch-and-rinse). DEI® New Extra Bond Two can also be used in a universal way with any type of composite made mainly of methacrylates, light-curing, self-curing and dual-curing.

Thanks to its specific formulation and in the presence of silica nanoparticles, this adhesive system of the latest generation

allows the creation of durable bond with the hard dental tissues by light-curing or dual curing (if previously mixed with the appropriate catalyst, DEI® New Extra Catalyst).

DEI® New Extra Bond Two is indicated for bond strength and longevity consistently high (> 30 MPA) dentin and enamel, regardless of the etching protocol used. The new functional monomer acid GDMA-P created for DEI® New Extra Bond Two guarantees a stable and lasting relationship with its chemical structure characterized by a hydrophobic chain, which ensures a good resistance to hydrolytic degradation, linked to functional groups that allows a strong phosphoric chemical bond with calcium ions present in the dental hydroxyapatite, without engaging the operator in complex sequences of applications.

This special formula allows DEI® New Extra Bond Two to be employed in steps cementing of indirect restorations (onlays, inlays, veneer), bridges and crowns through dual and self-curing resin systems, without firstly having to polymerize the adhesive layer. This is often a decisive factor when you want to reach difficult areas with the curing light.

WAY OF USE

- [A1] Adhesive for direct restorations in cavities of all classes with light-curing composite methacrylate-based.
- [A2] Adhesive for light-curing composite materials for core build-ups.
- [A3] Adhesive for indirect restorations with light-curing composites for cementation.
- [B] Intraoral repair of composite fillings, compomers, veneers and ceramic restorations or full zirconia without additional primers.
- [C] Treatment of exposed dentin with dentinal hypersensitivity.
- [D] Protective coating for glass ionomer cement fillings.
- [E1] Cavity sealing prior to amalgam restorations.
- [E2] Cavity sealing and core preparations before temporary cementation of indirect restorations.
- [F1] Adhesive for direct restorations in self-curing or dual-curing composite in all cavity classes and for core build-up.
- [F2] Adhesive for indirect restorations with self-curing or dual-curing composite for cementations.
- [G] Adhesive cementation of root posts with self-curing or dual-curing composite for cementations.

DEI® New Extra Bond Two creates a secure bond to both the natural dental hard tissue (enamel and dentin) and to based methacrylate composites, metals, silica-based ceramics and zirconia.

It is necessary to isolate the operative site from moisture areas. Avoid any contamination of the cavity or the surface of

the material with blood or saliva. We recommend the use of rubber dam.

The subsequent steps of the application depend on the type of clinical indication (see sections [A] - [G]).

In the case of light-curing materials, DEI® New Extra Bond Two can be applied without the dual catalyst DEI® New Extra Catalyst (indications [A] - [E]).

Mixing DEI® New Extra Bond Two + DEI® New Extra Catalyst is only required when using self-curing materials or with dual cure (Indications [F] - [G]).

[A] Direct or indirect restorations with light-curing composites for fillings, core build-ups and cementation.

1. **Prepare the cavity** or the core according to the principles of the adhesive technique used for fillings. Remove plaque and tartar in the event of non-carious Class V cavities. Rinse with plenty of water. Remove excess of water with a gentle stream of air. Do not dry the dentin. In case of deep cavities, protect the areas near of the pulp with a suitable sub-bottom of cavity based on calcium hydroxide.
2. **Etching of enamel** (and dentin, optional) with DEI® Etching Gel.
 - DEI® New Bond Two reaches high levels of adhesive strength on enamel and dentin when applied as a two-steps self-etching adhesive. Perform the etching of the dental hard tissue with DEI® Etching Gel (15-20 sec. in enamel) (5-10 sec. in dentin, optional).
 - It is still suggested a selective etching of the enamel at the beginning of the accession process in self-etching mode. In this case, apply DEI® Etching Gel on the prepared enamel and leave for 15-20 sec. Withdraw DEI® Etching Gel, rinse with water for about 10/15 sec. Remove the excess moisture with a gentle stream of air. Do not dry the dentin: the aim is to have a slightly moisture bright surface and not dry matte. Avoid excessively wet surfaces (over-wet).
3. **Apply DEI® New Extra Primer** thoroughly surfaces of the cavity and distribute it uniformly with a disposable applicator scrubbing for 15 sec. Dry with a jet of dry and oil-free air at maximum power on dental surfaces for 5 sec.
 - Repeat the application of DEI® New Extra Bond Two a second time and dry with a gentle stream of air on dental surfaces for 5 sec. in order to eliminate all the solvent.
4. DEI® New Extra Bond Two is a dual-curing adhesive and therefore can be used as dual-curing system immediately after mixing with catalyst (DEI® New Extra Catalyst). In this case, mix with an applicator in a mixing block, a drop DEI® New Extra Bond Two with a drop of DEI® New Extra Catalyst for about 3 sec. (the result is an etching dual-curing adhesive).

5. **Apply DEI® New Extra Bond Two carefully on the surfaces of cavities and distribute it evenly with a single-use applicator scrubbing for 5 sec. Dry with a stream of dry and oil-free air for 5 sec. in order to evenly distribute and remove excess adhesive from the cavity.**
6. Light cure the adhesive layer for 30 sec. with a device for polymerization (LED or halogen lamp with an emission > 500 mW/cm²).
7. Apply the restorative material for core build-up or cementation according to its instructions.

[B] Intraoral repair of fillings or restorations in composite, ceramic, metal (including amalgam).

1. Thoroughly roughen the surface of the restoration to which should be applied the adhesive and sandblast, preferably with aluminum oxide < 50 µm. In the case of ceramics, remove 1 mm of paint over the edge.
2. Clean thoroughly with air/water spray.
3. Continue with the application of DEI® New Extra Bond Two as previously described (see [A] 3. - 6.). It is not required the use of an additional primer.
4. Run the repair with a light-curing composite according to the manufacturer's instructions.

[C] Treatment of dentinal hypersensitivity on exposed dentin.

1. Clean the neck of the tooth or root surface with pumice powder and water or with a air-polishing system, preferably with glycine powder.
2. Clean thoroughly with water spray and dry with air free of moisture and oil. Do not over dry. The aim is to have a bright looking surface. Avoid over wet situations.
3. Next, continue with the application of DEI® New Extra Bond Two without using the etching with phosphoric acid (see [A] 3. - 6.).
4. Remove the layer of oxygen inhibition with a gauze or with cotton coiler.

[D] Protective coating for glass ionomer cement fillings.

1. If a matrix was used to facilitate the modelling, first of all it is necessary to remove it.
2. Apply a drop of DEI® New Extra Primer on a mixing block.
3. Apply DEI® New Extra Primer on the modelled filling, using a single use applicator for 5 sec.
4. Dry with dry and oil-free air for at least 5 sec. in order to remove the solvent.
5. Apply a drop of adhesive DEI® New Extra Bond Two on the modelled filling using a disposable applicator for 3 seconds.
6. Dry with dry and oil-free air for at least 5 sec. in order to

remove the excess adhesive.

7. Light cure DEI® New Extra Bond Two on every surface of filling for 20 sec. with a device for polymerization (LED or halogen lamp with an emission > 500 mW/cm²).
8. Once the filling has been completely polymerized, adjust the occlusion and seal once again areas instrumented with DEI® New Extra Bond Two in order to prevent the restoration from drying out.
9. Remove the layer of oxygen inhibition with a gauze or with a cotton roll.

[E] Cavity sealing before amalgam restorations or before the temporary cementation.

Note: do not seal the cavity with DEI® New Extra Bond Two if you want to make a temporary cementation using a material based on methacrylate.

1. Prepare the cavity according to the treatment plan. Clean with water. Wipe the excess water with a delicate air jet. Do not dry the dentin. The aim is to have a brilliant surface. Avoid over wet situations.
2. Protect the areas close to the pulp with a suitable liner cavity based on calcium hydroxide.
3. Apply the primer DEI® New Extra Primer and the adhesive DEI® New Extra Bond Two as explained in point [A] 3. - 8.
4. Light-cure on each surface of the restoration for 20 sec. with as described previously.
5. Continue with the application of the chosen material, according to the instructions of use.

[F] Direct or indirect restorations with self-curing or dual-curing composites for fillings, core build-ups and cementations.

The combination DEI® New Extra Bond Two + DEI® New Extra Catalyst creates a secure system for the application of self-curing and dual-curing composites, and methacrylate-based materials for cementations.

1. Prepare the cavity or the core according to the principles of adhesive technique for fillings. Clean with water. Dry the excess of water with a gentle stream of air. Do not dry the dentin. The aim is to have a bright looking surface. Protect areas near the pulp with a suitable liner of cavities based on calcium hydroxide.
2. Etching of the dental hard tissue (optional).
In self-etch mode, DEI® New Extra Bond Two reaches high levels of bond strength to enamel and dentin. If you want, you can run before the etching of dental hard tissue as previously explained in [A].
Etch enamel in cementing bridges Maryland. The etch self-curing adhesive materials are less effective on unprepared

enamel. The application of adhesive material on untreated or blunted enamel can cause discoloration or beveled margins. In case of restoration of Classes III, IV and V and diastema closures, for optimal results with DEI® New Extra Bond Two, enamel must be previously etched with DEI® Etching Gel.

- Etch surfaces with DEI® Etching Gel (according to the instructions of use).
 - Withdraw DEI® Etching Gel, rinse with water for about 20 sec. Remove excess moisture with a gentle stream air. Do not dry the dentin. The aim is to have a bright-looking surface. Avoid excessively moist surfaces (over-wet).
3. Apply DEI® New Extra Primer as explained before in [A].
 4. Thoroughly mix a drop of DEI® New Extra Bond Two and a drop of DEI® New Extra Catalyst on a mixing block with a disposable applicator for 3 sec.
 5. Apply the adhesive accurately on the surfaces of cavities and dispersing it for 20 sec. with an applicator.
 6. Dry the adhesive layer of DEI® New Extra Bond Two with dry and oil-free air for at least 5 sec. in order to remove the excess adhesive (See [A] for more information on the adhesive application).
 7. Light cure the adhesive layer for 20 sec. with a device for polymerization (LED or halogen lamp with an emission > 500 mW/cm²).
 8. According to the indication, continue with the application of chosen material, according to the instructions of use.

[G] Cementation of root canal posts.

1. Apply DEI® New Extra Primer for 20 sec. with an applicator in the root canal prepared before and on the occlusal surface as explained in [A].
2. Dry with dry and oil-free air for at least 5 sec. with maximum power. Repeat application and remove the solvent gently with a jet of air for 5 sec. (See [A] for more guidance on the application of the adhesive).
3. Thoroughly mix a drop of DEI® New Extra Bond Two and a drop of DEI® New Extra Catalyst on a mixing block with a single use applicator for 3 sec.
4. Remove excess material by blotting with a cone absorbent paper. Do not light cure the adhesive layer.
5. Apply adhesive for 5 sec. with an applicator in the root canal prepared before and on the occlusal surface.
6. Dry with dry and oil-free air for at least 3 sec. to remove excess adhesive from the canal.
7. Insert the root post in the canal, previously prepared according to the relevant instructions of the manufacturer, using a special self-curing or dual-curing composite for cementing the post (see the instructions for use) in order to obtain a

slight excess of material.

8. Cure DEI® New Extra Bond Two applied on occlusal, the material must be light cured for at least 40 seconds.

WARNINGS / PRECAUTIONS

- The phenolic substances (such as temporary materials that contain eugenol) prevent proper polymerization and therefore must not come into contact with DEI® New Extra Bond Two.
- Extra DEI® New Bond Two contains acidic monomers, UDMA, HEMA, camphorquinone, amines (DABE), BHT, catalysts and ethanol. Do not use DEI® New Extra Bond Two in case of hypersensitive (allergic) to any of these components.
- Avoid contact with the oral mucosa. Contact with DEI® New Extra Bond Two can cause bleaching of the oral mucosa as a result of coagulation of proteins. This is a reversible effect that resolves in a few days.
- In case of contact of DEI® New Extra Bond Two with eyes, flush eyes immediately with plenty of water and seek medical advice.
- Extra DEI® New Bond Two is light-cured. So avoid the excessive lighting during application.
- The mixed adhesive is self and light-curing and so it should be used immediately.

CONSERVATION

Close the cap of DEI® New Extra Primer, DEI® New Extra Bond Two and DEI® New Extra Catalyst immediately after using and when opened store at temperatures between 4°C-12°C.

Do not expose to direct sunlight or the operating light.

Do not use beyond the expiration date.

Our products have been created for use in the dental field. Regarding the application of our products, the verbal and/or written information is provided in accordance with our skills and with no obligation. Technical comments do not relieve the operator from carefully consult the documentation as to their suitability for the intended purposes. Since the use of DEI® italia products is not under our control, the operator is fully responsible for the clinical use. DEI® italia still guarantees the quality of the products distributed according to existing standards and in correspondence to conditions set out in the terms of sale and delivery.

Via Torino, 765 - 21020 Mercurio (Va) - Italy
Phone +39 0331 969270 - Fax +39 0331 969271
www.deitalia.it

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