## ZIRCOCAD 4D - DISCO FRESABILE IN ZIRCONIA

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## ZIRCOCAD 4D - ZIRCONIA MILLABLE BLOCKS

### Before use

ZircoCad 4D ESTHETIC and ZircoCad 4D INTEGRAL blocks, are pre-sintered dental zirconia. Please handle it carefully. Inspection each shipment for damage and do not use damaged discs for the production of dental restorations. Store ZircoCad 4D in a cool, dry, temperature stable environment (between 5°C and 50°C) in the original package.

CHARACTERISTICS	
ESTHETIC	INTEGRAL
Suitable for anterior and posterior restoration. Multilayer color gradient. High fracture toughness to resist chipping.	Suitable for anterior and posterior restoration. High strength, superior mechanical property. Color gradient, no need to do the dyeing process.

INDICATIONS	
ESTHETIC	INTEGRAL
Anterior restoration. Posterior restoration. Telescopic crown. 2 to 6 units bridge and more than 7 units bridge (Only one pontic).	Full contour crown and bridge. Maryland bridge. Telescopic crown.

TECHNICAL DATAS		
ZircoCad 4D	ESTHETIC	INTEGRAL
Sintering temperature	6.9 g/cm <sup>3</sup>	6.9 g/cm <sup>3</sup>
Flexural strength	700-1200 MPa	1100 MPa
Light transmittance	43-55%	42%
Sintering temperature	1500°C	1530°C

### RECOMMENDED DESIGN STEPS

#### SCANNING AND DESIGNING

Please scan with highly precise scanners to get accurate data of the restoration model. Design according to the conditions of patients and doctor's requirements.

#### MILLING

Pre-sintered zirconia material has an inherent shrinkage rate associated with each production.

The shrinkage rate must be input into the milling preparation software to ensure the accuracy of the eventual restoration.

When milling our zirconia, always follow these general guidelines:

- Please identify the incisal side at first.
- Only use sharp end mills with carbide or diamond coating.
- Do not use any restoration that has chips or cracks. Remove the units from the disc using a handpiece with a diamond-coated bur.
- Smooth the support areas with a medium-grit rubber polishing wheel.
- Remove any residual zirconia dust with an art brush.
- If a wet mill is used make sure all the zirconia is completely dry before sintering. Air dry for at least 15 minutes prior to sintering. Damp zirconia will crack if placed in the sintering oven.

#### SEPARATING RESTORATION AND CLEANING

Using technician specialized hand piece and grinder to separate restorations from the block.

Please clean up the powder on surface and inner side of restorations with brushes. If cleaning is not thorough, powder residues will contaminate coloring liquid when the dyeing is processing and the powder residues will stay on the surface and inner side of restorations after high temperature sintering, forming white spots and therefore having negative effect on esthetics of the restoration.

#### COLORING

Generally coloring is not needed for ZircoCad 4D multilayer zirconia.

#### SINTERING

Put the restoration in the crucible with zirconium beads. The bottom of crucible should be covered by the zirconium beads. Make sure the sintering furnace is clean.

Sinter the restoration according to the recommended temperature (see the sintering curve at the end attachment).

Recommended max heating rate: 10°C per minute.

Cooling should be done without temperature control in closed furnace.

Never open the furnace until temperature down to 80°C to avoid thermal shocks.


#### GRINDING AND POLISHING

Using special zirconia grinder to trim the surface of restoration. Steps: coarse grinding - fine grinding - rough polishing.


#### GLAZING

After coarsely finishing, proceed to external dyeing and glazing.

### ZIRCOCAD 4D ESTHETIC - SINTERING CURVE

### ZIRCOCAD 4D INTEGRAL - SINTERING CURVE

**EC REP** Prolinx GmbH - Brehmstr. 56,  
40239, Duesseldorf, Germany

**De Corematrix Co.Ltd** - No.40 Nanbeigang Road,  
Hukou County Jiujiang City Jiangxi China

**Changsha Honchon Technology Co. Ltd**  
International industry park, HN province, 410200, China

**DISTRIBUTED BY** Via Torino 765 - 21020 Merccallo (VA) Italy  
**DEI italia** Phone +39 0331 969270 - www.deiitalia.it

**CE**  
0197