DEI® CLEVER EASY-FILL

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DEI® CLEVER EASY-FILL

HIGHLY RESISTANT FLUID
COMPOSITE

Product for dental use only
Keep away from children.
Date of specifications: 03/2019

AVAILABLE PACKAGES

Two 2 g syringes & 10 tips
in the following colors:

A1 ^{REF} EFA1, A2 ^{REF} EFA2,
A3 ^{REF} EFA3, A3.5 ^{REF} EFA3.5,
A4 ^{REF} EFA4,
Traslucente ^{REF} EFTL,
White Opaque ^{REF} EFWO

EASY-FILL

CHARACTERISTICS

DEI® Clever Easy-Fill is a light curing fluid composite, radiopaque (210% Al), high viscosity composite. Complies with the guidelines described in ISO 4049. Its high percentage of filler offers better guarantees as to duration and resistance when used for permanent aesthetic restorations in areas undergone to loading.

COMPOSITION

- Diurethane dimethacrylate
- Butanediol dimethacrylate
- DOES NOT CONTAIN Bis-GMA

TECHNICAL DATA

Vickers hardness: 794 MPa
Flexural strength: 144 MPa
E-Module: 12316 MPa
Compressive strength: 431 MPa
Filler content: approx. 77%
Radiopacity: 210% Al

INDICATIONS

- Seal of fissure and extended fissure on molars and premolars.
- Fillings in class V cavities (cervical caries, eroded areas).
- Minimally invasive fillings in class I, II and III cavities.
- Restoring defects in enamel.
- Blocking out undercuts.
- Minimal adjustments to the contours and shade of the enamel.

ADVANTAGES

- Both flowable and viscous
- No running or dripping
- Easy and highly precise application.
- Even penetrates into areas that are difficult to reach.
- Stable and long-lasting restorations
- X-Ray check possible
- High polishability
- Addition in interproximal contacts
- Fill small bubbles or gaps
- Finish line (more elastic) on fiberglass or carbon fiber substrates
- Optimized biocompatibility

EASY-FILL WHITE OPAQUE

CHARACTERISTICS

DEI® Clever Easy-Fill White Opaque Baseline is a light curing fluid composite, highly radiopaque, does not contain Bis-GMA, suitable for hypersensitive patients. Marks the cavity bottom and its white opaque shade provides an optimal contrast to the filling material and natural tooth structure.

COMPOSITION

- Diurethane dimethacrylate
- Butanediol dimethacrylate
- DOES NOT CONTAIN Bis-GMA

TECHNICAL DATA

Vickers hardness: >654 MPa
Flexural strength: >114 MPa
E-Module: >10815 MPa
Compressive strength: >298 MPa
Filler content: approx. 73,1%

INDICATIONS

- Cavity bottom
- Defines the line between composite and tooth structure, and facilitates the work in case of retreatment
- Suitable for allergic patients because does not contain Bis-GMA
- Easy to find the same dark colour of natural tooth.

ADVANTAGES

- No dripping
- Easy and highly precise to apply
- Good wetting of tooth structure
- Maximum adhesion and seal
- For durable fillings
- High abrasion resistance
- Easy to polish
- Easy to detect on X-Ray
- Optimized biocompatibility
- Converts post-operative sensitivities and secondary caries
- Stain

WAY OF USE (STUDIO)

Preparatory measures
Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

1. Cavity preparation

Prepare the cavity as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Avoid preparations for finishing. Rinse the cavity with water, remove all debris and dry. The cavity must be isolated. It is advisable to use rubber dam.

2. Protection of the pulp

In very deep cavities, the areas closest to pulp must be coated with calcium oxide.

3. Approximal contact areas

For cavities with approximal sections, place the composite in a transparent matrix and fix it in place.

4. Etching

Apply etching gel absolutely precisely, beginning with the enamel margins. Allow it to react for at least 30 seconds before spraying thoroughly with water and drying the etched surfaces with oil-free compressed air. Avoid dehydrating the dentine. After drying, the etched surfaces of the enamel appear chalky/white and must not be contaminated prior to applying the bonding agent. Should they be contaminated with saliva, rinse and dry them again, re-etch if necessary.

5. Applying the bonding agent

Refer to the DEI® New Extra Bond instructions for details.

6. Applying DEI® Clever Easy-Fill

Apply the composite (max. 1,5 mm) of DEI® Clever Easy-Fill directly into the cavity using the curved application tips supplied with the material. When applying the material with the tip, ensure that no air bubbles become entrapped. Light cure each layer with a commercially available polymerization unit for 40 seconds, holding the light guide as close as possible to the surface and must not be touched or removed.

Finishing

Before curing the last layer of composite, apply DEI® Surface Hardener and photopolymerize for 40 seconds, wash an dry. The hardening of the composite shall be complete. The use of this product avoids the need of eliminating the layer of surface inhibition with burrs or other products, saving working time. If increase of aesthetics and restoration resistance is wanted, apply DEI® experience Seal Coat Fast. Photopolymerize for 60 seconds and leave cure for other 20 seconds.
Take out the surface layer of DEI® experience Seal Coat Fast using DEI® Polish Paste and small goat hair brush and polish

with Polish Felt.

At the end of the work an extremely hard surface, resistant and bright, is obtained. In addition, this technique considerably reduces the development of tartar and dental plaque.

NOTE

When placing time consuming restorations, to prevent premature curing of composite the light of the lamp should be temporarily moved away from the site temporarily of the composite covered with foil impervious to light.

For hygienic reasons, the curved application tips supplied with the material must only be used once.

Apply the composite material stratifying the various masses. Photo-polymerize every layer for 40 seconds for layers not thicker than about 1,5 mm, using photo-polymerizing lamps with an emission spectrum of at least 500 mW/cm². Should be impossible to bring sufficiently close the light source to the material, it is necessary lengthen the exposure time.

Halogen lamps are preferable as they offer wider light spectrum compared to LED lamps, so that guaranteed efficiency on any photo-sensitive material. To comply with curing times and layers is fundamental for achieving complete conversion of the materials, consequently resistance and duration in the time. The composites harden also after very few seconds, however they are not polymerized. In this phase saving time can cause a failure of the restoration.

WAY OF USE (LAB)

DEI® Clever Easy-Fill has chemical adhesion even on already light-cured composite.

Procedure:

Roughen the already light-cured composite, degrease with 96° ethyl alcohol, wet with DEI® Seal Coat Fast. Wait 1-2 minutes, do not light-cure, add DEI® Clever Easy-Fill and light cure.

STORAGE

Do not store above 25°C. Avoid direct sunlight. Close the syringe immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the continuous pressure let the product spill blocking the dispenser. Do not use after the expiry date.

This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

CONTRAINDICATIONS

Unpolymerized composite may cause skin allergies. The user must take adequate precautions (e.g. gloves). Should the operator or patient be known to be allergic to one of the constituents listed under "Composition", do not use this material.

SIDE-EFFECTS

To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide).

INTERACTION WITH OTHER SUBSTANCES

As phenolic substances (such as eugenol) inhibit polymerization, do not use products containing such substances.

Via Torino, 765 - 21020 Mercallo (VA) - Italy
Phone +39 0331 969270
Fax +39 0331 969271
www.deiitalia.it

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