DEI® RAINBOW EASYTEMP

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DEI® RAINBOW EASYTEMP

DUAL CURE COMPOSITE FOR TEMPORARY BRIDGES AND CROWNS

CE
0546

Product for dental use
Keep away from children
Date of the specifications: 02/2019

Available packages

One 75 g cartridge in the colours:

- A2 (REF. EASY 2)
- A3 (REF. EASY 10)
- 50 mixing cannulas (REF. EASY 4)
- Dispenser gun (REF. EASY 5)

Characteristics

DEI® Rainbow Easytemp is a self mixing, dual curing material consisting of two components. Shows fluorescence in UV-light. It is indicated for production of temporary crowns and bridges.

It offers very important advantages:

1. It swiftly and easily meets the requirements of aesthetic and functional temporary items.
2. The patient feels gratified by the prompt operation.
3. Working time is only few minutes (only 1-2 minute in the mouth).
4. After extraction from oral cavity it maintains during 5-6 minutes

- an elastic consistency useful for easier removal and elimination of eventual surplus.
5. As the duration of work is determined by self-mixing and capacity of photopolymerization, it is very easy to use.
6. It doesn't heat in the hardening phase (38°C).
7. Highly resistant to abrasion and flexure.
8. Top protection of the preparation.
9. Best polishing.
10. In case, repairable using the same material as a glue, after treatment of the two parts to connect with ethylic alcohol and DEI® Seal Coat Fast (not cured).

Composition

Polyfunctional acrylates and polyfunctional methacrylates; Unsaturated esters, malonyl urea derivatives barium glass silica filler.

Technical data

Resistance to compression: 270 MPa
Resistance to flexure: 70 MPa

Way of use

1. Take the impression either before the preparation of the stumps, or before the extraction; eliminate the interdental details from the impression.
2. Control that the outlets on the cartridge are free from hardened material; eventually, remove it using a probe. The material must be extruded slowly and continuously.
3. Insert the cannula on the cartridge rotating fully. The cannula shall re-

- main inserted on the cartridge up to the next use, acting as a top.
4. After discard of the first 2 mm, distribute the product inside the impression; place the impression inside the oral cavity within 30" and leave it there 1-2 minutes.
5. Extract the provisional that will be elastic for other 3-4 minutes and proceed to elimination of excesses and proximal undercuts. Wanting immediate hardening, photopolymerize every element for 60 sec with halogen light.
6. Take out the surface inhibiting layer with ethylic alcohol or acetone.
7. Finish using burs or rubber tops.
8. To get excellent aesthetic performance, apply DEI® Seal Coat Fast (not included in the package) and photopolymerize. Take out the surface part using abrasive paste and goat hair brush.
9. Fix with DEI® Easycem Hard.

Cautions

In case of individual sensitivity to product's components don't use.

Storage and shelf life

Don't use after the expiration date. Keep at temperature not beyond 20°C.

 **DEI® italia**
Via Torino, 765
21020 Mercallo (Va)
Phone +39.0331.969270
Fax +39.0331.969271
www.deiitalia.it

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